Dial Reading Viscometer

our original...over 75 years!

The World Wide Standard Viscometer

Easy to Select Speeds

Electronic Drive Means quiet, reliable operation

Analog Display

- Shows % Torque
- Use Factor Finder to convert reading to centipoise

Simple to use, easy setup

2 Year Limited Warranty

Available in explosion proof U.L. Class 1, Group D locations (w/o Electronic Drive)

Accuracy: ±1.0% of range

Repeatability: ±0.2%

Optional Accessories

Viscosity Standards

RV/HA/HB-1 Spindle

Quick Action Lab Stand

Temperature Bath

Small Sample Adapter

UL Adapter

Thermosel

Helipath Stand with T-bar Spindles

Spiral Adapter

DIN Adapter

Quick Connect/Extension Links

Vane Spindles



What's Included?

Instrument

6 spindles (RV/HA/HB) or 4 spindles (LV)

Spindle Guard Leg*

Lab Stand (Model A)

Carrying Case

*Not applicable to HA or HB models

VISCOSITY RANGE SPEEDS cP(mPa•s) MODEL Number o Min. Max. RPM LVT 1^{+} 2M .3-60 8 RVT 100† 8M .5-100 10 HAT 200++ 16M .5-100 10 HBT 64M 800†† .5-100 10

† 1 cP achieved with UL Adapter accessory. 15 cP on LV with standard spindles. †† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle. M = 1 million cP = Centipoise mPa•s = Millipascal•seconds



Easy Speed Adjustment and On/Off Control

